CEOS LAURIL SULFAT BUJON**NAMJENA**

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CEOS LAURYL SULFATE BROTH (LAURYL TRYPTOSE BROTH LTB)

USES

CEOS Lauryl Sulfate Broth (LTB) is liquid medium for the detection of coliforms in waters.

PRINCIPLES AND PROCEDURES

CEOS Lauryl Sulfate Broth (Lauryl Tryptose Broth – LTB) is a selective medium recommended for the enumeration of coliforms in water and dairy products as well as for confirmatory tests of lactose fermentation with gas production by coliforms in foods. Another advantage of this medium is that the indole test can be performed directly in the tube. APHA recommends the use of Lauryl Tryptose Broth for the Most Probable Number Presumptive Test of coliforms in waters, effluent or sewage, also as a confirmatory test of lactose fermentation with gas production for milk samples, and for the detection of coliforms in foods. This broth was elaborated to promote a rich growth and high gas production from small inocula of coliform organisms. Lauryl Sulfate Broth is recommended by the norm ISO 4831 and ISO 7251 for the detection and enumeration of coliforms and E. coli respectively, by the most probable number technique. The coliform group is both aerobic and anaerobic. Tryptose provides the nitrogen and vitamins required for organism growth. Lactose is the fermentable carbohydrate. Sodium chloride maintains the osmotic balance of the medium. Sodium lauryl sulphate is the selective agent used to inhibit organism other than coliforms. Potassium phosphates are the buffering agents.

REAGENTS

CEOS Lauryl Sulfate Broth (g/l)

Sodium chloride	5.00
Lactose	5.00
Tryptose	20.00
Potassium dihydrogen phosphate	2.75
Sodium lauryl sulfate	0.1
Dipotassium hydrogen phosphate	2.75

pH 6.8±0.2

Color: light amber.

PRECAUTIONS

IVD. Only for professional use. Do not use product if there are visible signs of microbial contamination, color change, drying, cracking or other signs of quality deterioration.

Apply aseptic techniques and identified precautions against microbiological hazards throughout all procedures.

Prepared plates, bottles, sample bottles and other contaminated materials have to be sterilized in autoclave after use and before disposing.

STORAGE AND SHELF LIFE

Products that are listed on cardboard packaging are packed in tight plastic foil, for better stability. Products must be kept at dark place and in temperature span suggested on labels, in original packaging until use. Keeping at temperature below 5°C and repeated and/or extreme differences in temperatures can amplify moisture formation. Moisture in form of fine mist or small drops on the inside of the lid, especially during and after keeping in the fridge, is acceptable and sign of medium freshness. Products can be used until expiration date (see packaging label) and incubated during recommended incubation periods.

If products from opened stack are kept on clean place where temperature is between 2-8°C, they can be used for one week (7 days). Exposure to light before and during incubation must be reduced.

Shelf life: 6 months.

USERS QUALITY CONTROL

Inoculate representative samples with following type cultures:

Microorganismis	Results
<i>Enterococcus faecalis</i> ATCC 19433	No growth
<i>E.coli</i> ATCC 25922	Gas production and turbidity
<i>Enterococcus faecalis</i> ATCC 29212	No growth
<i>Citrobacter freundii</i> ATCC 43864	Gas production and turbidity
<i>E.coli</i> ATCC 8739	Gas production and turbidity

According to ISO 11133:

- For detection of coliforms ISO 4831
Incubation conditions: Productivity, Selectivity: 24±2 h to 48±2 h / 30±1 °C.
Inoculation conditions: Productivity qualitative: <100 CFU, Selectivity: 10⁴-10⁶ CFU
- For detection of E.coli ISO 7251
Incubation conditions: Productivity, Selectivity: 24±2 h to 48±2 h / 37±1 °C.
Inoculation conditions: Productivity qualitative: <100 CFU, Selectivity: 10⁴-10⁶ CFU

PROCEDURE

Materials provided

CEOS Lauryl Sulfate Broth. Microbiologically controled.

Materials not provided

Plates for extra cultures, reagents and laboratory equipment if needed.

Test procedure

For the detection and enumeration of coliforms according to ISO 4831:

- Inoculate the tubes of selective enrichment with the test portion.
- For inoculum of 1 ml or less, use single strength LTB medium.
- For inoculum of 10 ml or more, use double strength LTB medium.
- Incubate at 30 °C or 37 °C (as agreed) for 24 h (double strength medium) or, for 24 h and another additional 24 h if neither gas formation nor opacity preventing the detection of gas formation is not observed at this stage (single strength medium).
- From the incubated tubes of double strength LTB and incubated tubes showing gas formation or opacity of single strength LTB, inoculate a tube of confirmation medium (Brilliant Green Bile Broth) and observe after 24 or 48 h.
- A tube with gas formation is a positive tube.
- Count the total number of positive tubes.

For the detection and enumeration of presumptive *E. coli* according to ISO 7251:

- Inoculate the tubes of selective enrichment broth (Lauryl Tryptose Broth) with the initial suspension.
- For inoculum of 1 ml or less, use single strength LTB medium.
- For inoculum of 10 ml or more, use double strength LTB medium.
- Incubate the tubes of LTB at 37 °C for up to 48 h, and examine the gas production after 24 h and 48 h.

- Each tube of double-strength LTB that has given rise to opacity, cloudiness or gaseous emission, and each tube of single-strength LTB that has given rise to gaseous emission, is subcultured to a tube containing EC Broth.
- Incubate the tubes of EC Broth at 44 °C for up to 48 h, and examine the gas production after 24 h and 48 h.
- Each tube of EC Broth that has given rise to gaseous emission is subcultured to a tube containing Indole-Free Peptone Water and incubated at 44° C for 48 h.
- Tubes showing opacity, cloudiness or gas production in LTB and whose subcultures have produced gas in EC broth and indole in Peptone water, are considered as positive tubes that containing presumptive E coli.
- Count the total E. coli by the MPN method.

PACKAGING/AVAILABILITY

Tubes, ready to use, 20x10 ml	531310
-------------------------------	--------

MORE INFORMATION

For more information contact manufacturer.

C-PHARM d.o.o.

Gospodarska zona 15

32000 Vukovar

Tel.: +385-32-534-154

Mail: info@c-pharm.hr

Web: www.c-pharm.hr